



STONE BASKET

Helical and Flat Wire Stone Basket

Helical Wire



Flat Wire



Issuing Body: COM - May '04 Last Revision: '00

GTA STONE BASKET retrieval is particularly designed for treatment of ureter stones.

Available also in Nitinol wire.

Nitinol wire offers a kink resistant, flexible design for maximum scope deflection.

Stainless Steel

| Order Code Flat Wire | Order Code Helical Wire | Description | Unit per Box |
|----------------------|-------------------------|----------------|--------------|
| JB 0370 | TV 0370 | 3 Fr-70 cm-4W | 1 |
| JB 0390 | TV 0390 | 3 Fr-90 cm-4W | 1 |
| JB 0312 | TV 0312 | 3 Fr-120 cm-4W | 1 |
| JB 0470 | TV 0470 | 4 Fr-70 cm-4W | 1 |
| JB 0490 | TV 0490 | 4 Fr-90 cm-4W | 1 |
| JB 0412 | TV 0412 | 4 Fr-120 cm-4W | 1 |
| JB 0570 | TV 0570 | 5 Fr-70 cm-4W | 1 |
| JB 0590 | TV 0590 | 5 Fr-90 cm-4W | 1 |
| JB 0512 | TV 0512 | 5 Fr-120 cm-4W | 1 |

Nitinol

| Order Code Flat Wire | Order Code Helical Wire | Description | Unit per Box |
|----------------------|-------------------------|----------------|--------------|
| ROM 0370 | MIL 0370 | 3 Fr-70 cm-4W | 1 |
| ROM 0390 | MIL 0390 | 3 Fr-90 cm-4W | 1 |
| ROM 0312 | MIL 0312 | 3 Fr-120 cm-4W | 1 |
| ROM 0470 | MIL 0470 | 4 Fr-70 cm-4W | 1 |
| ROM 0490 | MIL 0490 | 4 Fr-90 cm-4W | 1 |
| ROM 0412 | MIL 0412 | 4 Fr-120 cm-4W | 1 |
| ROM 0570 | MIL 0570 | 5 Fr-70 cm-4W | 1 |
| ROM 0590 | MIL 0590 | 5 Fr-90 cm-4W | 1 |
| ROM 0512 | MIL 0512 | 5 Fr-120 cm-4W | 1 |

Features:

- Available in both versions, helical and flat wire
- Wire available in 3, 4, 5 French
- Polyethylene sheath
- On request it is possible to produce needles of different size