

# PICA

## Pigtail Universal Drainage Catheter

### PICA: Drainage Set with Catheter & Needle



### PICA L: Drainage Set with Stringlock Catheter & Needle



### Pigtail Tip



Issuing Body: COM - May '04 Last Revision: '00

### PICA: Drainage Set with Catheter & Needle

Code	Code catheter only	Size (F x cm)	Units per box
10010520	30010520	5 x 20	5
10010530	30010530	5 x 30	5
10010620	30010620	6 x 20	5
10010630	30010630	6 x 30	5
10010720	30010720	7 x 20	5
10010730	30010730	7 x 30	5
10010820	30010820	8 x 20	5
10010830	30010830	8 x 30	5
10010930	30010930	9 x 30	5
10011030	30011030	10 x 30	5
10011230	30011230	12 x 30	5

### Features:

- Radiopaque polyurethane catheter
- Stiffened cannula with trocar stylet
- Needle assembly locks to catheter during insertion
- Choice of 20 cm or 30 cm catheter from 5 F to 8 F
- Wide choice of sizes from 5 F up to 12 F (6 F - 12 F Stringlock)
- Pica L: Pigtail catheter with stringlock for added security
- Catheter only option for over the wire placement
- Large oval holes in pigtail for maximum drainage

### PICA L: Drainage Set with Stringlock Catheter & Needle

Code	Code catheter only	Size (F x cm)	Units per box
15010620	35010620	6 x 20	5
15010630	35010630	6 x 30	5
15010720	35010720	7 x 20	5
15010730	35010730	7 x 30	5
15010820	35010820	8 x 20	5
15010830	35010830	8 x 30	5
15010930	35010930	9 x 30	5
15011030	35011030	10 x 30	5
15011230	35011230	12 x 30	5

# PICA-HYDROPHILIC

## Hydrophilic Coated Pigtail Drainage Catheter

### PICA H: Drainage Set with Catheter & Needle



### PICA HL: Drainage Set with Stringlock Catheter & Needle



### Pigtail Tip



Issuing Body: COM - May '04 Last Revision: '00

### Features:

- Radiopaque polyurethane catheter, hydrophilic coated
- Stiffened cannula with trocar stylet
- Needle assembly locks to catheter during insertion
- Choice of 20 cm or 30 cm catheter from 5 F to 8 F
- Wide choice of sizes from 5 F up to 12 F (6 F - 12 F Stringlock)
- Pica HL: Pigtail catheter with stringlock for added security
- Catheter only option for over the wire placement
- Large oval holes in pigtail for maximum drainage

### PICA H: Drainage Set with Catheter & Needle

Code	Code catheter only	Size (F x cm)	Units per box
10110520	30110520	5 x 20	5
10110530	30110530	5 x 30	5
10110620	30110620	6 x 20	5
10110630	30110630	6 x 30	5
10110720	30110720	7 x 20	5
10110730	30110730	7 x 30	5
10110820	30110820	8 x 20	5
10110830	30110830	8 x 30	5
10110930	30110930	9 x 30	5
10111030	30111030	10 x 30	5
10111230	30111230	12 x 30	5

### PICA HL: Drainage Set with Stringlock Catheter & Needle


Code	Code catheter only	Size (F x cm)	Units per box
15110620	35110620	6 x 20	5
15110630	35110630	6 x 30	5
15110720	35110720	7 x 20	5
15110730	35110730	7 x 30	5
15110820	35110820	8 x 20	5
15110830	35110830	8 x 30	5
15110930	35110930	9 x 30	5
15111030	35111030	10 x 30	5
15111230	35111230	12 x 30	5

# URONEPHRO

## Pigtail Nephrostomy Set - Double Puncture

### The set include:

- Two-piece nephrostomy needle with the possibility of kidney puncture with the 22G needle and the coaxial introducing of the 18G needle;
- PTFE coated flexible guide wire with a 3mm J tip;
- Four dilators length 20 cm;
- 30 cm drainage connecting tube;
- Pigtail catheter, with 1-way female Luer Lock stopcock;
- Retention disk;
- Disposable scalpel;
- 2000 cc drainage bag.



### ORDERING INFORMATION:

Order code	French	Length (cm)	G.W. (inch)	# pcs. per box
31070530	5	30	.028"	1
31070640	6	40	.035"	1
31070740	7	40	.035"	1
31070840	8	40	.038"	1
31070940	9	40	.038"	1
31071040	10	40	.038"	1

### Only for replacement catheter

Order code	# pcs. per box
31010530	10
31010640	10
31010740	10
31010840	10
31010940	10
31011040	10

**GTA**<sup>®</sup>  
INTERNATIONAL  
MEDICAL DEVICES


GTA srl  
via della Vittoria, 5  
46026 Quistello (MN) Italy  
Tel. +39 0376 625300  
Fax +39 0376 625279  
E-mail: gta@mynet.it

# URONEPHRO

## Reduced Pigtail Nephrostomy Set - Double Puncture

### The set include:

- Two-piece nephrostomy needle with the possibility of kidney puncture with the 22G needle and the coaxial introducing of the 18G needle;
- PTFE coated flexible guide wire with a 3mm J tip;
- Dilator in the same french-size as the drainage catheter; the catheter is made with a pigtail curvature;
- 30 cm drainage connecting tube;
- Pigtail catheter, with 1-way female Luer Lock stopcock;
- Retention disk.



### ORDERING INFORMATION:

Order code	French	Length (cm)	G.W. (inch)	# pcs. per box
31050530	5	30	.028"	5
31050640	6	40	.035"	5
31050740	7	40	.035"	5
31050840	8	40	.038"	5
31050940	9	40	.038"	5
31051040	10	40	.038"	5

### Only for replacement catheter

Order code	# pcs. per box
31010530	10
31010640	10
31010740	10
31010840	10
31010940	10
31011040	10

**GTA**<sup>®</sup>  
INTERNATIONAL  
MEDICAL DEVICES

GTA srl  
via della Vittoria, 5  
46026 Quistello (MN) Italy  
Tel. +39 0376 625300  
Fax +39 0376 625279  
E-mail: gta@mynet.it