


BOHR

Bone Marrow Biopsy/Aspiration Needle



Issuing Body: COM - May '04 Last Revision: '00

Features:

- Ergonomic handle design
- Cannula with 5 cutting edges to facilitate penetration during biopsy
- Stylet with trocar point for a smoother entry into the marrow
- Available with the explant model with 5 lateral holes on the external cannula and an additional stylet

BOHR Standard version

Code BOHR	Size (gaxcm)	Units per Box
MN 0810	08 x 10	5
MN 0815	08 x 15	5
MN 0910	09 x 10	5
MN 0915	09 x 15	5
MN 1110	11 x 10	5
MN 1115	11 x 15	5
MN 1307	13 x 07	5
MN 1310	13 x 10	5

BOHR Special version for Explant

Code BOHR Explant Model	Size (gaxcm)	Units per Box
MN 1310	13 x 10	5

On request it is possible to produce needles with different sizes

NAXI

Disposable Sternal Iliac Aspiration Needle



Issuing Body: COM - May '04 Last Revision: '00

Special Model for Explant



2 lateral holes for an increased aspiration

Features:

- Adjustable and removable needle stop to control depth of penetration
- Luer lock connector
- Millimetric feed
- Colour code for easy gauge identification

NAXI Standard version

Code NAXI	Size (gaxmm)	Units per Box
NAXI 1501	15 x 10-44	10
NAXI 1503	15 x 25-60	10
NAXI 1601	16 x 10-44	10
NAXI 1603	16 x 25-60	10
NAXI 1801	18 x 10-44	10
NAXI 1803	18 x 25-60	10

NAXI Special version for Explant

Code NAXI Explant Model	Size (gaxmm)	Units per Box
NAXI 1404	14 x 37-71	10

On request it is possible to produce needles with different sizes