


UROSTENT

Ureteral Double Pigtail Stent

Ureteral Double Pigtail Stent



Issuing Body: COM - May '04 Last Revision: '00

Complete Set (including G.W. and clamp) - Open Tip

Length \ Fr	16 cm	22 cm	24 cm	26 cm	28 cm	30 cm
4.7	20034716	20034722	20034724	20034726	20034728	20034730
5	20030516	20030522	20030524	20030526	20030528	20030530
6	20030616	20030622	20030624	20030626	20030628	20030630
7	20030716	20030722	20030724	20030726	20030728	20030730

Complete set (including G.W. and clamp) - Closed Tip

Length \ Fr	16 cm	22 cm	24 cm	26 cm	28 cm	30 cm
4.7	20044716	20044722	20044724	20044726	20044728	20044730
5	20040516	20040522	20040524	20040526	20040528	20040530
6	20040616	20040622	20040624	20040626	20040628	20040630
7	20040716	20040722	20040724	20040726	20040728	20040730

Single Catheter - Open Tip

Length \ Fr	16 cm	22 cm	24 cm	26 cm	28 cm	30 cm
4.7	20134716	20134722	20134724	20134726	20134728	20134730
5	20130516	20130522	20130524	20130526	20130528	20130530
6	20130616	20130622	20130624	20130626	20130628	20130630
7	20130716	20130722	20130724	20130726	20130728	20130730

Single Catheter - Closed Tip

Length \ Fr	16 cm	22 cm	24 cm	26 cm	28 cm	30 cm
4.7	20144716	20144722	20144724	20144726	20144728	20144730
5	20140516	20140522	20140524	20140526	20140528	20140530
6	20140616	20140622	20140624	20140626	20140628	20140630
7	20140716	20140722	20140724	20140726	20140728	20140730


Features:

- Available in two basic configurations: double open pigtail and double closed pigtail
- Wide choice of sizes from 4.7 F up to 7 F, and from 16 cm to 30 cm catheter
- The complete set includes:
 - PTFE coated guide wire, 100 cm long and fixed core
 - Pusher
 - polyurethane stent with side holes along the entire catheter
 - Clamp
- Catheter only option for over the wire placement

UROSOFTSTENT

Ureteral Double Pigtail Stent


Ureteral Double Pigtail Stent



Guide wire



Pusher



Clamp



Issuing Body: COM - May '04 Last Revision: '00

Complete Set (including G.W. and clamp) - Open Tip

Length \ Fr	16 cm	22 cm	24 cm	26 cm	28 cm	30 cm
4.7	22034716	22034722	22034724	22034726	22034728	22034730
5	22030516	22030522	22030524	22030526	22030528	22030530
6	22030616	22030622	22030624	22030626	22030628	22030630
7	22030716	22030722	22030724	22030726	22030728	22030730

Complete set (including G.W. and clamp) - Closed Tip

Length \ Fr	16 cm	22 cm	24 cm	26 cm	28 cm	30 cm
4.7	22044716	22044722	22044724	22044726	22044728	22044730
5	22040516	22040522	22040524	22040526	22040528	22040530
6	22040616	22040622	22040624	22040626	22040628	22040630
7	22040716	22040722	22040724	22040726	22040728	22040730

Single Catheter - Open Tip

Length \ Fr	16 cm	22 cm	24 cm	26 cm	28 cm	30 cm
4.7	22134716	22134722	22134724	22134726	22134728	22134730
5	22130516	22130522	22130524	22130526	22130528	22130530
6	22130616	22130622	22130624	22130626	22130628	22130630
7	22130716	22130722	22130724	22130726	22130728	22130730

Single Catheter - Closed Tip

Length \ Fr	16 cm	22 cm	24 cm	26 cm	28 cm	30 cm
4.7	22144716	22144722	22144724	22144726	22144728	22144730
5	22140516	22140522	22140524	22140526	22140528	22140530
6	22140616	22140622	22140624	22140626	22140628	22140630
7	22140716	22140722	22140724	22140726	22140728	22140730

Features:

- Available in two basic configurations: double open pigtail and double closed pigtail
- Wide choice of sizes from 4.7 F up to 7 F, and from 16 cm to 30 cm catheter
- The complete set includes:
 - PTFE coated guide wire, 100 cm long and fixed core
 - Pusher
 - Soft polyurethane stent with side holes along the entire catheter
 - Clamp
- Catheter only option for over the wire placement