

# GALAK

## Galactography Procedure Kit



### This Kit includes:

- 1 cc syringe;
- 25 gauge blunt needle (curved) with a transparent Luer Lock;
- Set of blunt round dilators: 23 gauge-yellow, 24 gauge-dark green, 26 gauge-light green, 30 gauge-orange, and 31 gauge-light blue.

### ORDERING INFORMATION:

Order code	# pcs. per box
------------	----------------

GALAK 10



# GALA

## Galactography Set

### The Set includes:

- 27 gauge blunt cannula with transparent Luer Lock
- 15 cm polyethylene tube with outer diameter 1.25 mm and inner diameter 0.75mm
- Round dilator with outer diameter of 0.25mm and 0.30mm (ordering code GALAD)


### ORDERING INFORMATION:

Order code	# pcs. per box
------------	----------------

GALA 10

GALAD 10

Sterile Single Use Sets



**IGTA**<sup>®</sup>  
INTERNATIONAL  
MEDICAL DEVICES

GTA srl  
via della Vittoria, 5  
46026 Quistello (MN) Italy  
Tel. +39 0376 625300  
Fax +39 0376 625279  
E-mail: gta@mynet.it

# SCIA

## Sialography Set

This product is used for the cannulation of salivary ducts, submandibular ducts, and parotid ducts.


### The Set includes:

- 27 gauge blunt needle with transparent Luer Lock;
- 2mm (outer diameter) PVC tube with an inner diameter of 1 mm

### ORDERING INFORMATION:

Order code	# pcs. per box
SCIA	10

Sterile Single Use Needle



# LYM

## Lymphography Set

LYM is used for the injection of contrast media into the lymphatic vessels.


### The Set includes:

- 1,5cm sharpened cannula with transparent Luer Lock;
- 60cm polyethylene tube.

### ORDERING INFORMATION:

Order code	Outer d. x inner d. of tube (mm)	Gauge of cannula (cm)	# pcs. per box
LYM 27	2.00 x 1.00	27	10
LYM 30	0.70 x 0.30	30	10

Sterile Single Use Needle



**IGTA**<sup>®</sup>  
INTERNATIONAL  
MEDICAL DEVICES

GTA srl  
via della Vittoria, 5  
46026 Quistello (MN) Italy  
Tel. +39 0376 625300  
Fax +39 0376 625279  
E-mail: gta@mynet.it