


ANGIOPLASTY BALLOON CATHETER

Angioplasty PTA Balloon Catheter - .035i and .018i

Guide Wire 0.18i					
Balloon Diameter mm	Balloon Length mm	Introducer F	Rated Burst Pressure	Catalogue Number Usable Length 90 cm	Catalogue Number Usable Length 120 cm
2,0	20	4	20	-	94202020
	40	4	20	-	94202040
2,5	20	4	20	-	94225020
	40	4	20	-	94225040
3,0	20	4	20	-	94203020
	40	4	20	-	94203040
	100	4	18	-	94203010
4,0	20	4	20	94104020	94204020
	40	4	20	94104040	-
	60	4	18	94104060	-
5,0	20	4	20	94105020	-
	40	4	20	94105040	-
	80	4	16	94105080	-
6,0	20	4	17	94106020	-
	40	4	17	94106040	-
	60	4	14	94106060	-
7,0	20	4	14	94107020	-
	40	4	14	94107040	-
	60	4	14	94107060	-
8,0	20	5	14	94108020	-
	40	5	14	94108040	-
9,0	20	5	13	94109020	-
	40	5	13	94109040	-



Issuing Body: COM - May '04 Last Revision: '00

OUTSTANDING TRACKABILITY

- Short tapered tip design
- Highly kink-resistant shaft design
- Low deflated profile;
- Excellent pushability

EXCELLENT DILATION RESULTS

- Uniform pressure distribution;
- Highly scratch and puncture resistant
- High rated burst pressure

Guide Wire 0.35i					
Balloon Diameter mm	Balloon Length mm	Introducer F	Rated Burst Pressure	Catalogue Number Usable Length 75 cm	Catalogue Number Usable Length 120 cm
4,0	20	5	20	93104020	93204020
	40	5	20	93104040	93204040
	60	5	20	93104060	93204060
	100	5	20	93104010	93204010
5,0	20	5	18	93105020	93205020
	40	5	20	93105040	93205040
	60	5	20	93105060	93205060
	80	5	18	93105080	93205080
6,0	20	5	18	93106020	93206020
	40	5	19	93106040	93206040
	60	5	19	93106060	93206060
	80	6	18	93106080	93206080
7,0	20	5	17	93107020	93207020
	40	5	16	93107040	93207040
	60	6	16	93107060	93207060
	80	6	17	93107080	93207080
8,0	20	5	17	93108020	93208020
	40	5	15	93108040	93208040
	60	6	15	93108060	93208060
9,0	20	6	15	93109020	93209020
	40	6	15	93109040	93209040
	60	6	15	93109060	93209060
10	40	6	15	93110040	93210040
	60	7	13	93110060	93210060
12,0	20	7	14	93112020	-
	40	7	12	93112040	-
14,0	20	7	10	93114020	-
	40	7	10	93114040	-