


# **HYSKA**

## **Hysterosalpingography Catheter**



HYSKA facilitates injection of contrast media into the uterine cavity to evaluate patency of the fallopian tubes.

The catheter design allows good physician control while maintaining sufficient flexibility for comfort.

Issuing Body: COM - December '04 - Last Revision: '03

### **Features:**

- The rounded catheter tip and large side holes assures safe placement and superior flow
- Choice of two sizes 5 F and 7 F
- Balloon latex free
- The catheter is supplied with introducer, stopcock and syringe

Code	Size (Fr x cm)	Units per Box
HSC 530	5 x 30	5
HSC 730	7 x 30	5